

# CATALOG 2015

BRIEF EXPLANATION OF THE  
IMPROVEMENTS AND NEW  
MODELS

# FIXED HEIGHT TREATMENT COUCHES

- CD-F3

3 sections and 3-level adjustable fixed height treatment couch with detachable legs.



- CM-20  
Wood couch with backrest

# 2 SECTIONS COUCHES

- CE-2129-AR

2 section electric examination couch,  
retractable wheels, breathing hole plug



## CE-2129-R

2 section electric examination couch with  
negative adjustment of head section,  
retractable wheels, breathing hole plug



### Characteristics:

Structure with rounded shapes

Head section with gas piston

Maximum weight of 180 kg

Important: the standard price include hand switch and  
don't include breathing hole plug

# 2 SECTIONS COUCHES

- CE-0128 (SPRING PREVIEW)  
Especially recommended for aesthetic

## Characteristics:

X shape structure

Electric height elevation via hand switch

Vertical elevation (without annoying horizontal movement)

Mechanical adjustment of head section

Head section with breathing hole



# 3 SECTIONS COUCHES

## ■ CE-2139-R

3 section electric examination couch with retractable wheels and negative adjustment of head section.

Structure with rounded shapes

Head and foot section with gas piston



## • CE-0138

3 section electric examination couch.

Structure with rounded shapes

Head and foot section with mechanical adjustment.



Important: the standard price include hand switch and don't include birthing hole plug

# 9 SECTION COUCHE

- CE-2399-R
- 9 section multifunctional couch with motorised elevation, backrest, Trendelenburg tilt and drainage positioning adjustments, with retractable weels.

Electric height elevation via hand switch  
Adjustment by gas piston of head section  
Motorised trendelemburg tilt 30°.  
Motorised drainage facility between mid and backrest section  
Flat fold-down armrests  
Head section with breating hole and breating hole plug  
Height adjustable backrest armrests  
Pedal-based retractable wheel design



# IMPROVEMENTS

Rounded shape of the table



Modification of structure of models ended en 7  
(e.g. CE/CH 0127).

The shape of pedal for rise the whells

The shape of the base of the structure is more esthetical

The anchorage of the whells